

# Morita Scans Viewer Manager

Just put the CD into the CD drive and click here to import a new CBCT scan into the manager.

Use this if the scan data is not given to you on a CD

Select your patient's name

Click here to open the scan in 3D

Just type in your comments relative to each image. Scroll through your previous notes too.

A reminder of the scanning parameters used for each scan

Patient ID	Patient Name
00000006489	
0000000693	

Photo Date	Open XYZ Viewer	Open 3D Viewer
11/05/2011		

Date	Image Comments
12/05/2011	
12/05/2011	

Image info

kV:90  
mA:5.5  
PRJAGL:1PI  
EXPTIME:17.4  
PIXEL:80.0um  
VOLUME\_RADIUS:20.000  
RECON\_FILTER:D\_001

Export the selected images to the patient's notes

Click to make thumbnail larger

At a later date, review your screenshots as a reminder of the specificities of the case, without having to go back to the full 3D data



CAVENDISHIMAGING